**Bürkle GmbH**

Rheinauenstr. 5 | 79415 Bad Bellingen

Telefon +49 7635 6121-0

info@buerkle.de

www.buerkle.de

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## General Safety Instructions

- ▶ Read these operating instructions carefully before use and keep them accessible for all users.
- ▶ The device may only be used by trained and qualified personnel.
- ▶ Observe the manufacturer's specifications as well as the safety and hazard instructions of the medium used.
- ▶ Use the device only within the limits of material compatibility and for its intended purpose.

## Safety Instructions Before Use

- ▶ Use only with approved liquids that are compatible with the materials in contact with the product.
- ▶ Do not use with aggressive or easily flammable/combustible liquids without separate approval.
- ▶ Wear appropriate personal protective equipment, at minimum safety goggles and protective gloves. Use additional protective clothing depending on the medium.
- ▶ For indoor use only. Place the source and receiving containers securely and ensure the hose cannot slip out unintentionally.
- ▶ Before starting the transfer, check that the receiving container has sufficient capacity.
- ▶ Before each use, inspect the siphon for visible damage, leaks, secure connections, and proper function of the ball valve.
- ▶ Protect the device from prolonged UV exposure.

## Function

The siphon is designed for sampling and transferring liquids from higher containers to lower containers according to the siphon principle.

It is available for sale in two versions:

- ▶ 5605-1000 with flexible tube
- ▶ 5605-1100 with rigid tube

## Operating Limits

- ▶ Permissible operating temperature: 5 °C to 35 °C.
- ▶ Permissible viscosity: max. 500 mPas.
- ▶ If the medium or application changes, suitability must be reassessed.
- ▶ Do not continue using the device if it is damaged or malfunctioning.

## Scope of delivery

1 x	Siphon with stopcock	①
1 x	Directions for use	②

You can find additional accessories at [www.buerkle.de](http://www.buerkle.de)

## Operation

1. Put on personal protective equipment.
2. Immerse the long suction hose end into the source container until it is in the liquid. Route the hose without kinks and secure it against slipping out.
3. Close the ball valve. Do not place fingers in the movement area of the ball valve, as there is a slight risk of pinching.
4. Repeatedly squeeze and release the suction bulb until liquid is drawn in.
5. Open the ball valve. The liquid will flow automatically into the lower container.
6. Monitor the fill level. Do not leave the transfer process unattended.
7. To stop, close the ball valve, completely empty the hose, and remove the device with minimal dripping.

### CAUTION!

- ▶ After priming, the liquid flow starts automatically.
- ▶ Liquid may escape due to overflow, dripping losses, leakage, or a slipping hose.
- ▶ Risk of slipping, contamination, and injury due to the medium used.
- ▶ Monitor the fill level, change the receiving container in time, and properly close the ball valve after filling

## Technical Data

- ▶ Flow rate: approx. 9 l/min
- ▶ The device consists of a suction hose and discharge hose, hanging bend (PP), hose nozzles (PP), suction bulb (LDPE) and union nuts (PP), as well as a ball valve (PP + NBR) and flat gaskets (LDPE). Depending on the version, the suction hose and discharge hose are made of PVC (5605-1000) or LDPE (5605-1100).

## After Use / Maintenance

- ▶ Rinse the device with a suitable cleaning medium after use.
- ▶ Clean thoroughly when changing media to prevent cross-contamination and product carryover.
- ▶ Store the siphon clean, dry, and without kinks, and protect it from UV exposure.
- ▶ If used frequently, regularly check for leaks, hose aging, and proper function of the ball valve.
- ▶ Replace damaged parts immediately

## Troubleshooting

Problem	Solution
Suction bulb does not draw in liquid	Close the ball valve. Check the suction hose for blockages. Ensure the hose is sufficiently immersed in the liquid.
Liquid does not flow	Check whether the ball valve is open and the suction bulb is sufficiently filled with liquid. Repeat the priming process if necessary.
Flow rate too low	Fully open the ball valve. Check hose routing and remove kinks. Observe permissible viscosity.
Liquid leaks	Check connections, hose, and ball valve for tightness and damage. Do not continue using the device if damaged.
Receiving container overflows	Supervise the transfer process, monitor the fill level, and close the ball valve in time.
Hose slips out of the container	Position the hose correctly and secure it against slipping out unintentionally.
Dripping after transfer	Fully close the ball valve and empty the hose after the transfer process.
Cross-contamination / residual medium	Clean the device thoroughly and rinse completely when changing media.

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