

**0309-0xxx** Pressure sprayer, various size, food

## Declaration of Conformity for food contact material

**Our supplier hereby confirms that the above mentioned plastic product complies with the Regulations, Directives and Laws as described below.**

### 1. General

This certificate has been issued in accordance with Regulation (EC) No. 1935/2004 of the European Parliament and of the Council on materials and articles intended to come into contact with food.

With regard to the migration tests carried out, the samples tested comply with the following legislation:

Article 3 (1) a), b) of Regulation (EC) No. 1935/2004 of the European Parliament and of the Council on materials and articles intended to come into contact with food Commission Regulation (EU) No. 10/2011 on plastic materials and articles intended to come into contact with food. Sections 30 and 31 (1), (2) a) and (3) of the German "Food, Commodities and Feed Code" (LFGB) in conjunction with the "Commodities Ordinance" (BedGgstV).

### SML-substances

2. Substances used for which restrictions and/or specifications are in place under Regulation (EC) No. 10/2011 incl. amendments

- For the investigated product all contained substances with a specific migration limit (SML) were compiled (based on the datasheets relevant for food-contact of all used raw materials) and the compliance with the limit values ascertained by analyses, calculations, or any other plausible way. - FCM 19, 20, 132, 234, 282, 368, 433, 779

### Non-intentionally added substances (NIAS)

- Potentially present NIAS (impurities, degradation- and reaction-products), were investigated via a GC/MS-screening. In case NIAS were detectable, consumer safety was ascertained by internationally accepted principles for risk-assessment.

### Inertness of the food contact material

- A confirmation that the organoleptic properties (smell, taste and colour) of foods are not changed unacceptably was part of the test scope to prove compliance with Article 3 of the Framework-Regulation (EC) No 1935/2004 as well as the investigation of overall migration (if legally required for the tested article).

### **Dual-Use-additives**

3. Substances which are subject to a restriction in food (Dual-use-additives) - E171, E210, E470a, E471, E553b

### **Specifications regarding the use of the plastic, in particular:**

4.1 Intended type(s) of food(s)

- All types of food

# ZERTIFIKAT CERTIFICATE

Pumpen | Probennehmer | Laborbedarf  
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware  
for Laboratory, Industry, Science



**0309-0130** Pressure sprayer, 1200 ml, food

4.2 Intended conditions for food contact

- Any long-term storage at room-temperature or below incl. heating to (2 hrs/ 70 °C or 15 min/ 100 °C)

4.3 Tested surface/Volume-ratio

- ONE ARTICLE / 1 dm<sup>2</sup> with 180 ml food simulant

**Use of a functional barrier accordin to 10/2011 Article 3**

- none

The printing ink is not food-compatible, but the printing on the pressure pump sprayer is harmless, as it is mounted on the outside of the atomizer and therefore it does not come into contact with the contents.

Further tests are not carried out by us.

The customer must test whether the product is suitable for the application.

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of certificate is 002284-QM15.

September 4, 2026

A handwritten signature in blue ink, appearing to read "M. Saint-Denis".

Bürkle GmbH, Bad Bellingen  
Martin Saint-Denis, Managing Director