

0310-xxxx Wash bottle, LDPE transparent, various size

Declaration of Compliance

Food Contact Compliance Assessment

Compliance with General Food Contact Legislation

Bürkle GmbH confirms that the products listed above are manufactured from low-density polyethylene (LDPE) and, based on the available compliance documentation, comply with the requirements of Regulation (EC) No 1935/2004 and Commission Regulation (EU) No 10/2011, including the amendments covered by the available supporting documentation, under the specified conditions of use.

The products are manufactured in accordance with the applicable principles of Good Manufacturing Practice pursuant to Commission Regulation (EC) No 2023/2006.

OML - overall migration limit

According to the available compliance documentation, the overall migration of the products complies with the applicable legal limits established by Commission Regulation (EU) No 10/2011 under the specified conditions of use.

SML - specific migration limit

According to the available compliance documentation, the applicable specific migration limits (SML) are complied with under the specified conditions of use.

Where substances subject to specific restrictions are present, the corresponding legal limits are met.

For substances not included in the Union List pursuant to Annex I of Commission Regulation (EU) No 10/2011, additional national provisions, such as the recommendations of the German Federal Institute for Risk Assessment (BfR), may need to be considered.

Food Contact Suitability

Based on the available compliance assessment, the products are suitable for long-term contact with all types of food at room temperature or below.

The following temperature and time conditions are also covered:

- Heating up to 70 °C for a maximum of 2 hours
- Heating up to 100 °C for a maximum of 15 minutes

The compliance assessment is based on a food contact surface area-to-food mass ratio of:

6 dm²/kg

No functional barriers are used.

ZERTIFIKAT CERTIFICATE

Pumpen | Probennehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



0310-xxxx Wash bottle, LDPE transparent, various size

No additional application-specific testing is performed by Bürkle GmbH. The user is responsible for verifying the suitability of the product for the intended application and conditions of use.

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of certificate is 002284-QM15.

October 21, 2026

A handwritten signature in blue ink, appearing to read "M. Saint-Denis".

Bürkle GmbH, Bad Bellingen
Martin Saint-Denis, Managing Director