

ZERTIFIKAT CERTIFICATE

Pumpen | Probennehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



0355-xxxx

Wide-necked box, square, PVC transp., various size, w/cap

Type(s) of food intended to come into contact with the material	All types of food groups in terms of the Regulation (EU) No. 10/2011
Type(s) of food not to come into contact with the material	None
Duration and temperature of treatment and storage while in contact with the food	Long-term-storage at up to 40°C
Ratio of the area of the surface coming into contact with the food to volume used to determine compliance of the specific total material or article	9,6 cm ² to 20 ml simulant
Specification of tested conditions of use with maximum test requirements (without exceeding OML)	3% Acetic Acid: 10 days 40°C 10% Ethanol: 10 days 40°C Olive Oil: 10 days 40°C
<p>The overall migration values are below the statutory limit (under the above-mentioned conditions of use). Testing is carried out in accordance with Articles 17 and 18 of Commission Regulation (EU) No. 10/2011 on plastic materials and articles intended to come into contact with food, in conjunction with Annex V. Testing at 40 °C for 10 days covers any storage duration under refrigerated and frozen conditions, including hot filling and/or heating to temperatures between 70 °C and 100 °C (70 °C < T < 100 °C) for a maximum of 2 hours at 70 °C and 15 minutes at 100 °C.</p>	
Information about restricted substances and / or specifications	<ul style="list-style-type: none"> • 2-Ethyl-1-hexanol: 30 mg/kg • Dodecyl acrylate: 0.05 mg/kg • Tin compounds: 0.18 mg/kg (expressed as tin, Sn) • Vinyl chloride: < 1 mg/kg • Butadiene: < 1 mg/kg
Information on Dual-Use Additives in accordance with Regulation (EU) No. 10/2011	Sodium salts of fatty acids (E470a)
Functional barrier	None
Further confirmation to the raw material	<p>BSE/TSE – Note for Guidance EMEA/410/01, Rev. 3 RoHS – Directive 2011/65/EU, including Directive (EU) 2015/863 REACH – Regulation (EC) No. 1907/2006, Annex II SVHC – Regulation (EC) No. 1907/2006, Annex XIV Phthalates – Not used in the manufacture or formulation of the product. Packaging and Packaging Waste – Directive 94/62/EC of 20 December 1994 on packaging and packaging waste, as amended by Directive 2005/20/EC and Directive 2013/2/EU Nanomaterials Swiss VOC Allergens</p>

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft
Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



0355-xxxx

Wide-necked box, square, PVC transp., various size, w/cap

NIAS	<p>A NIAS screening was performed using 95% ethanol as the food simulant.</p> <p>Substances listed in Annex I of Commission Regulation (EU) No. 10/2011 were not detected at concentrations exceeding the applicable limits. For any non-listed substances that may be present, a risk assessment has been carried out by the raw material manufacturer in accordance with Article 19 of Regulation (EU) No. 10/2011.</p>
------	---

The user must convince himself of the suitability of the product beyond the requirements of the guidelines.

Further tests are not carried out by us.

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of certificate is 002284-QM15.

July 15, 2026