

0533-0010

"Spigot ""Air"" self-venting, PE, 3/4"" outer- NW 13mm"

Food Packaging compliance Statement

We can confirm compliance to general legislations:

1. Regulation (EC) No 1935/2004 on materials and articles intended to come into contact with food.
2. Commission Regulation (EU) No 10/2011 on plastic materials and articles intended to come into contact with food and amendments.
3. Commission Regulation (EC) No 2023/2006 on good manufacturing practice for materials and articles intended to come into contact with food.

The silicone component (washer) is formulated to comply with the following requirements:

1. BfR Recommendation XV. For Silicones.
2. FDA CFR 21 § 177.2600 Rubber articles intended for repeated use.

Overall Migration Testing (for components in contact with food)

- The overall migration limit is less than 10mg/dm² (60mg/kg).
- The chosen test conditions are according to the requirements of EU Regulation No. 10/2011.
- Migration test: 100cm² of the sample / 100 ml of food simulant
- Food simulants:
 - Simulant A - 10% ethanol
 - Simulant B - 3% acetic acid
 - Simulant D2 - olive oil or alternate fatty food simulant (95% ethanol)
- Test conditions: 40°C / 10 days

Specific Migration Testing (for components in contact with food)

- The specific migration values of the substances and additives used are in compliance with migration limits required by EU Regulation No. 10/2011 (Annex I and Annex II).
- The chosen test conditions are according to the requirements of EU Regulation No. 10/2011.
- Migration test: 100cm² of the sample / 100 ml of food simulant
- Food simulants:
 - Simulant A - 10% ethanol
 - Simulant B - 3% acetic acid
 - Simulant D2 - olive oil or alternate fatty food simulant (95% ethanol)
- Test conditions: 60°C / 10 days

The product is considered to be suitable for food contact according to Regulation (EC) No. 1935/2004 of the European Parliament and of the Council. The articles shall not unacceptably change the organoleptic characteristics of food.

Type(s) of food intended to come into contact with the product: All types of food.

Intended food contact condition: Any long-term storage at room temperature or below, including heating up to 70°C for up to 2 hours, or heating up to 100°C for up to 15 minutes.

ZERTIFIKAT CERTIFICATE

Pumpen | Probennehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



0533-0010

"Spigot ""Air"" self-venting, PE, 3/4"" outer- NW 13mm"

The user must convince himself of the suitability of the product beyond the requirements of the guidelines. Further tests are not carried out by us.

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of certificate is 002284-QM15.

April 17, 2025

A handwritten signature in blue ink, appearing to read "M. Saint-Denis".

Bürkle GmbH, Bad Bellingen
Martin Saint-Denis, Managing Director