

7102-xxxx Graduated beaker w/handle, PP, various sizes, transp.grad.

Declaration of conformity

for plastic materials that come into contact with food.

Our supplier hereby declares that the products comply with the statutory provisions of the Consumer Goods Ordinance and the European Regulations (EC) No. 1935/2004 and (EC) No.10/2011 in their respective current versions.

Testing of the products by an independent, accredited institute for compliance with the limit values did not reveal any exceedances. The sensory tests also did not reveal any odor or taste impairments and were carried out in accordance with Art. 17 and 18 of Regulation (EU) No. 10/2011 in conjunction with Annex V. The materials and raw materials used therefore comply with the Consumer Goods Ordinance and Regulation (EC)

No. 10/2011 and may be used for the manufacture of food contact materials and articles with the restrictions specified therein with regard to the migration limits and permissible residual contents in the end product.

Furthermore, it can be confirmed that the raw material currently used to manufacture the above-mentioned products was produced according to the specifications of the raw material manufacturer in compliance with the requirements of good manufacturing practice (GMP) for products that may come into contact with food (in accordance with European Regulation (EC) No. 2023/2006).

If the products contain substances with restrictions (SML /QM), the limit values listed in Regulation (EC) No. 10/2011 and recommendations of the BfR are complied with.

The texts of the guidelines and recommendations can be downloaded in the current version from the Internet at <http://eur-lex.europa.eu> and <http://bfr.bund.de> respectively.

Specification for the intended use or restrictions:

- Types of food that may come into contact with the material:

All types of food (dry, aqueous, acidic, alcoholic, fatty) - tested according to Table 3 Annex III with the food simulants listed there (1. dist. water / water of the same quality or food simulant A (ethanol 10% vol. %); 2. food simulant B (acetic acid 3% by volume) and 3. food simulant D2 (any vegetable oil with less than 1% unsaponifiable matter).

- alternatively with 95 % ethanol and isooctane according to chapter 2 section 2.1.3 Contact conditions when using food simulants).

- Types of food that must NOT come into contact with the material:

(none)

- Duration and temperature of handling / storage in contact with the food:

Tested according to Table 3 Annex IV for all types of food, for 10 days at 40°C (isooctane for 2 days at 20°C) or for 10 days at 60°C - covers any storage period under refrigeration and freezing conditions, including hot-filling and/or heating to a temperature T, of $70^{\circ}\text{C} < T < 100^{\circ}\text{C}$, for a maximum duration of $t = 120/2^{((T-70)/10)}$ minutes.

- There are no restrictions on the starting materials used to manufacture the products

(SML values) are available:

SML: -mglkg

- Note on "dual-use substances":

The material does not contain any substances, according to Regulation (EC) 1333/2008 in its current version, which are also permitted as food additives.

No functional barrier made of plastic is used.

Traceability of the product according to Regulation (EC) No. 1935/2004 is ensured by the use of a date stamp on the product by the use of a date stamp on the product.

The user must satisfy himself of the suitability of the product for the intended the intended filling material.

The declaration is based on our current state of knowledge and the information provided by our supplier at the time the document was created.

The manufacturer, Bürkle GmbH in Bad Bellingen, has been certified accordance with the DIN EN ISO 9001 standard by DQS (German Association for Quality Assurance (DQS) since 1995, Certificate No. 002284-QM15.

20. March 2025



Bürkle GmbH, Bad Bellingen,
Martin Saint-Denis, Managing Director

