

7503-xxxx Grad.Cylinder PP blue sc. various size

Food Contact Statement

The products listed above comply with the applicable requirements for plastic materials and articles intended to come into contact with food.

The conformity assessment is based on testing and technical documentation covering the products specified in this certificate. The products fulfil the relevant requirements of:

- Regulation (EC) No. 1935/2004 on materials and articles intended to come into contact with food
- Commission Regulation (EU) No. 10/2011 on plastic materials and articles intended to come into contact with food (including applicable amendments)
- Commission Regulation (EC) No. 2023/2006 on good manufacturing practice (GMP)

Testing carried out by an independent accredited laboratory confirmed compliance with the applicable overall migration limits. Organoleptic testing showed no adverse influence on the taste or odour of food under the specified test conditions

Intended use

The products are suitable for contact with all food types, including:

- dry food
- aqueous food
- acidic food
- alcoholic food
- fatty food

The declaration applies under the test conditions specified in Commission Regulation (EU) No. 10/2011.

Test Results

Organoleptic Properties

Testing according to DIN EN 10955 confirmed compliance with the applicable requirements. No unacceptable changes in taste or odour were detected.

Overall Migration

Overall migration testing according to DIN EN 1186 demonstrated compliance with the maximum permitted migration limit of 10 mg/dm² for all specified food simulants

Overall migration testing was carried out according to DIN EN 1186 using the food simulants specified in Regulation (EU) No. 10/2011 under the prescribed test conditions.

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



7503-xxxx Grad.Cylinder PP blue sc. various size

Simulant	Unit	Measurement uncertainty	Sample	1. Contact	2. Contact	3. Contact = Measured value	Assessment
Acetic Acid 3 %	mg/dm ²	10 %	1.	< 1	< 1	< 1	passed
			2.	< 1	< 1	< 1	passed
			3.	< 1	< 1	< 1	passed
Ethanol 10 %	mg/dm ²	10 %	1.	< 1	< 1	< 1	passed
			2.	< 1	< 1	< 1	passed
			3.	< 1	< 1	< 1	passed
Olive oil	mg/dm ²	30 %	1.	< 4	< 4	< 4	passed
			2.	< 4	< 4	< 4	passed
			3.	< 4	< 4	< 4	passed
Isooctane	mg/dm ²	30 %	1.	3.3	3.2	3.1	passed
			2.	2.1	2.1	2.1	passed
			3.	2.9	2.7	2.6	passed
Ethanol 95 %	mg/dm ²	30 %	1.	< 1	< 1	< 1	passed
			2.	< 1	< 1	< 1	passed
			3.	< 1	< 1	< 1	passed

Specific Migration

Specific migration testing for the regulated substances confirmed compliance with the applicable limits.

Additional Information

The products do not contain functional barriers. Traceability is ensured by product marking or lot identification.

The declaration applies exclusively to the products listed above and under the intended conditions of use. Users remain responsible for verifying suitability for their specific application.

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



7503-xxxx Grad.Cylinder PP blue sc. various size

The user must convince himself of the suitability of the product beyond the requirements of the guidelines.

Further tests are not carried out by us.

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of certificate is 002284-QM15

July 16, 2026

A handwritten signature in blue ink, appearing to read "M. Saint-Denis".

Bürkle GmbH, Bad Bellingen
Martin Saint-Denis, Managing Director