

ZERTIFIKAT CERTIFICATE

Pumpen | Probennehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



8300-1010 Hand scoops for foodstuffs

8300-501x Filling scoops for foodstuffs

8300-001x 1-part shovels for foodstuffs

8300-2010 2-part shovel for foodstuffs

8300-401x Scrapers for foodstuffs

We hereby declare that the above products, material polymer: polypropylene, comply with the following regulations:

- Verordnung (EC) Nr 1935/2004
- Verordnung (EC) Nr 2023/2006
- Verordnung (EU) Nr 10/2011

The material has been tested in accordance with the requirements of Commission Regulation (EU) No 10/2011 on plastic materials and articles intended to come into contact with foodstuffs using the BSEN 1186:2002 method.

The regulations stipulate that no plastic must be capable of releasing its constituents into foodstuffs with which it may come into contact in quantities exceeding the corresponding limit value. For the material the appropriate limit is 10 mg/dm².

The migration from the material was less than the maximum permitted by the Regulations.

Good manufacturing practice:

All procedures regarding the manufacturing of these products, including raw material supply, storage, processing, quality control, testing and packing are in accordance, adhere to and are compliant with European Directive EU 2023/2006. In respect of European Commission regulation# 2023/2006 of 22 December 2006 on good manufacturing practice for materials and articles intended to come into contact with food. This particular regulation refers specifically to EC regulation 1935/2004/EC in terms of materials.

This is to confirm that the ingredients used to manufacture the products listed below and the way these materials are handled, the processes they are put through are all subject to the Quality Assurance system (ISO9001 :2008) as approved by ISOQAR. As such this means we, and the products listed below meet the European Commission regulation # 2023/2006.

This is to confirm that this master batch is formulated and manufactured using materials of a synthetic origin using good manufacturing practices that meet European Commission regulation # 2023/2006. There are no ingredients in the formulation of our hygiene PP material that is of animal origin. As such, this material should not pass on any Animal derived disease like BSE (Bovine Spongiform Encephalopathy) or other TSE (Transmissible Spongiform Encephalopathy).

We have to inform you that our hygiene material contains traces (1-10 ppm) of a phthalate, originated from the used catalyst system. These traces fully comply with the EC Directives 2005/84/EC and Commission Regulation (EU) 10/2011 and amendments.

We can also inform you, that this material is not subject to Annex XIV (Authorisation) of Regulation (EC) No 1907/2006 (also known as REACH), since the possible traces phthalate present in our material are either regarded as an impurity or are far below the threshold of 0.1% (1000 ppm) as mentioned in Article 56(6) (b) of REACH (see also our REACH declaration).

ZERTIFIKAT CERTIFICATE

Pumpen | Probennehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



Specifications for use:

Type or types of food with which it may come into contact: All types of food

Type of intended use of the product: Repeated, intermittent use

Duration and temperature of treatment and storage in contact with food:

Any long-term handling at room temperature or below, including up to 60°C for up to 2 hours.

Maximum short-term operating temperature between -30°C and +80°C.

Optimum long-term operating temperature between +5°C and +40°C.

Ratio between the surface area in contact with food and the volume to determine the conformity of the product: 2 dm²/1 dl.

All migration tests were carried out by a UKAS accredited testing laboratory.

Overall migration test results:

Olive oil, 4 hours at 20°C: <3 mg/dm²

95% ethanol, 24 hours at 40°C: <1 mg/dm²

Iso-octane, 4 hours at 20°C: <3 mg/dm²

3% acetic acid, 24 hours at 20°C: <2 mg/dm²

The overall migration limit is 10 mg/dm². The migration from the material was less than the maximum permitted by the Regulations. Specific migration of metals: Specific migration testing was carried out for barium, cobalt, copper, iron, lithium, manganese and zinc. All reported results complied with the applicable permissible limits. Overall result: PASS.

We confirm that no allergens that are recognised by the EU FIC and FDA are used in the production of our products.

Bisphenol A (BPA): The material of these products is not intentionally manufactured or formulated with 4,4'-Bisphenol A (BPA), CAS No. 80-05-7, EC No. 201-245-8 and is compliant with (EU) 2018/2013.

Polyfluoroalkyl substances (PFAS): These substances as such are not intentionally used or added.

Specific migration limits (SML) and dual-use additives:

According to the information provided by our supplier, the material contains no monomers regulated with a specific migration limit (SML) and no intentionally incorporated additives regulated with a specific migration limit (SML). The material contains one or more intentionally added dual-use additives.

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of certificate is 002284-QM15.

September 3, 2026

A handwritten signature in blue ink, appearing to read "M. Saint-Denis", written over a light blue horizontal line.

Bürkle GmbH, Bad Bellingen

Martin Saint-Denis, Managing Director