

**8700-xxxx** Connector straight, conical nozzle, PP, various size

## Declaration of Conformity

Food Contact Materials

### 1. General Compliance

Bürkle GmbH confirms that the above-mentioned polypropylene (PP) articles comply with the applicable requirements for materials and articles intended to come into contact with food.

The articles comply with the applicable requirements of Regulation (EC) No 1935/2004 and Commission Regulation (EU) No 10/2011, including subsequent amendments. According to the available conformity documentation, the articles have been subjected to migration testing and fulfil the applicable requirements of Commission Regulation (EU) No 10/2011.

### 2. Intended Food Contact and Conditions of Use

The conformity assessment covers repeated contact with all types of food. The articles are suitable for food contact under conditions including heating up to 70 °C for a maximum of 2 hours or up to 100 °C for a maximum of 15 minutes, provided that such heating is not followed by prolonged storage at room temperature or under refrigerated conditions.

### 3. Migration Testing

Migration testing was performed using the standardized test condition OM3 (2 hours at 70 °C) in accordance with Annex V, Chapter 3.1 of Commission Regulation (EU) No 10/2011. For the assessment of food contact suitability, a surface area-to-food ratio of 6 dm<sup>2</sup>/kg was applied.

### 4. Overall and Specific Migration

For the assessment of overall migration into fatty foods, appropriate alternative food simulants, including 95 % ethanol and isooctane, were used in accordance with Annex V, Chapter 3, Section 3.4.2 of Commission Regulation (EU) No 10/2011. The highest migration value was taken into account for the assessment. Specific migration testing was also performed using the appropriate alternative simulant.

### 5. Basis and Scope of this Statement

This statement is based on the documentation available to Bürkle GmbH and our current state of knowledge. It applies to the above-mentioned articles when used under the specified conditions. In the event of deviating application conditions, the user is responsible for verifying the conformity and suitability for the intended application.

# ZERTIFIKAT CERTIFICATE

Pumpen | Probennehmer | Laborbedarf  
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware  
for Laboratory, Industry, Science



**8700-xxxx** Connector straight, conical nozzle, PP, various size

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of certificate is 002284-QM15

October 5, 2026

A handwritten signature in blue ink, appearing to read "M. Saint-Denis".

Bürkle GmbH, Bad Bellingen  
Martin Saint-Denis, Managing Director