

Material certificate

The media-contacting parts of the articles on the following pages are manufactured from stainless steel or aluminum, some in combination with PTFE. They are all designed to have a smooth finish and to be internally free of internal crevices.

The declaration for stainless steel is based on REGULATION (EC) No 1935/2004 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL: "The principle underlying this regulation is that any material or article intended to come into contact, directly or indirectly, with food must be sufficiently inert to preclude substances from being transferred to the food in quantities large enough to endanger human health or to bring about an unacceptable change in the composition of the food or a deterioration in its organoleptic properties."

Stainless steel is widely used and generally suitable for most foodstuffs. If using highly corrosive or aggressive materials, please check our list of "Chemical-Resistance":

www.buerkle.de/en/chemical-resistance or contact us for suitability.

This declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up.

Further tests are not carried out by us.

The declaration is based on our current state of knowledge and the information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the DIN EN ISO 9001 standard by the DQS (German Society for Quality Assurance) since 1995. The certificate number is 002284-QM15.

June 22, 2026

A handwritten signature in blue ink, appearing to read "M. Saint-Denis".

Bürkle GmbH, Bad Bellingen
Martin Saint-Denis, Managing Director

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



Article	Description	Alu	PTFE	PA	1.4034	1.4301 304	1.4305	1.4310	1.4571 316Ti	1.4404 316L
0600-3100	Thread ad. AISI304 inner/inner G 1 1/2" - G 2"		x			x				
2309-0006	Bucket, stainless steel, AISI 304, capacity 6 l					x				
2309-0010	Bucket, stainless steel, AISI 304, capacity 10 l					x				
2309-0013	Bucket, stainless steel, AISI 304, capacity 13 l					x				
2309-0013	Bucket, stainless steel, AISI 304, capacity 13 l					x				
2309-0015	Bucket, stainless steel, AISI 304, capacity 15 l					x				
2601-1001	Safety can 316L, dosage spout, relief valve, 1 l		x						x	
2601-1002	Safety can 316L, dosage spout, relief valve, 2 l		x						x	
2601-1005	Safety can 316L, dosage spout, relief valve, 5 l		x						x	
2601-2010	Safety can 316L, dosage spout, relief valve, 10 l		x						x	
2601-2020	Safety can 304, dosage spout, relief valve, 20 l		x			x				
2601-3005	Safety jerrycan 316L, screw cap, relief valve, 5 l		x						x	
2601-3010	Safety jerrycan 316L, screw cap, relief valve, 10l		x						x	
2601-3020	Safety jerrycan 304, screw cap, relief valve, 20 l		x			x				
2601-4010	Safety storage cont., 316L, tab, relief valve, 10l		x						x	
2601-4025	Safety storage cont., 304, tab, relief valve, 25 l		x			x				
2601-5050	Safety storage container 304, relief valve, 50 l		x			x				
2601-5125	Safety storage container, 304, relief valve, 25 l		x			x				
2601-6010	Safety storage cont., 316L, tab, level indic., 10l		x						x	
2601-7010	Safety jerrycan 316L, w/o relief valve, 10 l		x						x	
4208-0001	Tray stainless steel, AISI 304, 280x180x20 mm					x				
4208-0002	Tray stainless steel, AISI 304, 330x230x20 mm					x				
4208-0003	Tray stainless steel, AISI 304, 460x300x20 mm					x				
4208-0004	Tray stainless steel, AISI 304, 520x395x20 mm					x				
4208-0005	Tray stainless steel, AISI 304, 300x200x20 mm					x				
4208-0006	Tray stainless steel, AISI 304, 450x350x20 mm					x				
4208-0018	Bowl stainless steel, shallow form, 150x100x40 mm					x				
4208-0019	Bowl stainless steel, shallow form, 200x100x40 mm					x				
4208-0020	Bowl stainless steel, shallow form, 200x150x50 mm					x				
4208-0021	Bowl stainless steel, shallow form, 280x180x50 mm					x				
4208-0022	Bowl stainless steel, shallow form, 300x200x50 mm					x				
4208-0023	Bowl stainless steel, shallow form, 300x250x50 mm					x				
4208-0024	Bowl stainless steel, shallow form, 400x300x50 mm					x				
4208-0025	Bowl stainless steel, shallow form, 500x300x50 mm					x				
4208-0026	Bowl stainless steel, shallow form, 230x120x50 mm					x				
4208-0028	Bowl stainless steel, shallow form, 250x200x50 mm					x				
4208-0031	Bowl stainless steel, shallow form, 400x200x50 mm					x				

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



Article	Description	Alu	PTFE	PA	1.4034	1.4301 304	1.4305	1.4310	1.4571 316Ti	1.4404 316L
4208-0032	Bowl stainless steel, shallow form, 450x300x50 mm					x				
4208-0040	Bowl stainless steel, deep form, 400x200x100 mm					x				
4208-0041	Bowl stainless steel, deep form, 400x300x100mm					x				
4208-0042	Bowl stainless steel, deep form, 500x300x100mm					x				
4709-0004	Mixing bowl/laboratory bowl, round, AISI304, 18 ml					x				
4709-0008	Mixing bowl/laboratory bowl, round, AISI304, 150ml					x				
4709-0010	Mixing bowl/laboratory bowl, round, AISI304, 300ml					x				
4709-0015	Mixing bowl/laboratory bowl, round, 304, 1000 ml					x				
4709-0020	Mixing bowl/laboratory bowl, round, 304, 2500 ml					x				
4709-0025	Mixing bowl/laboratory bowl, round, 304, 4500 ml					x				
4709-1004	Laboratory bowl, conical/flat bottom, 304, 125 ml					x				
4709-1005	Laboratory bowl, conical/flat bottom, 304, 330 ml					x				
4709-1007	Laboratory bowl, conical/flat bottom, 304, 530 ml					x				
4709-1010	Laboratory bowl, conical/flat bottom, 304, 940 ml					x				
4709-1012	Laboratory bowl, conical/flat bottom, 304, 1720 ml					x				
4709-1014	Laboratory bowl, conical/flat bottom, 304, 6650 ml					x				
5301-0500	Extension AISI304, 100 cm, SiloDrill					x				
5301-1000	SiloDrill, alu, 150 cm, spiral screw Ø 90mm, 400ml	x	x			x		x		
5301-2000	SiloDrill, 304, 150 cm, spiral screw Ø 90mm, 400ml		x			x		x		
5302-0301	PowderProof f.granules, lance horiz.30cm, complete							x		x
5302-0302	PowderProof f.powders, lance horiz. 30cm, complete							x		x
5302-0601	PowderProof granules, lance vert. 60 cm, complete							x		x
5302-0602	PowderProof powders, lance vertical 60cm, complete							x		x
5302-0901	PowderProof granules, lance vert. 90 cm, complete							x		x
5302-0902	PowderProof powders, lance vertical 90cm, complete							x		x
5302-1201	PowderProof granules, lance vert. 120 cm, complete							x		x
5302-1202	PowderProof powders, lance vertic. 120cm, complete							x		x
5307-0101	Head, volume 0,1 ml, MicroSampler tube-Ø 12 mm									x
5307-0102	Head, volume 0,2 ml, MicroSampler tube-Ø 12 mm									x
5307-0103	Head, volume 0,3 ml, MicroSampler tube-Ø 12 mm									x
5307-0104	Head, volume 0,4 ml, MicroSampler tube-Ø 12 mm									x
5307-0105	Head, volume 0,5 ml, MicroSampler tube-Ø 12 mm									x
5307-0106	Head, volume 0,6 ml, MicroSampler tube-Ø 12 mm									x
5307-0107	Head, volume 0,7 ml, MicroSampler tube-Ø 12 mm									x
5307-0108	Head, volume 0,8 ml, MicroSampler tube-Ø 12 mm									x
5307-0109	Head, volume 0,9 ml, MicroSampler tube-Ø 12 mm									x
5307-0110	Head, volume 1 ml, MicroSampler tube-Ø 12 mm									x

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



Article	Description	Alu	PTFE	PA	1.4034	1.4301 304	1.4305	1.4310	1.4571 316Ti	1.4404 316L
5307-0111	Head, volume 1,1 ml, MicroSampler tube-Ø 12 mm									x
5307-0112	Head, volume 1,2 ml, MicroSampler tube-Ø 12 mm									x
5307-0113	Head, volume 1,3 ml, MicroSampler tube-Ø 12 mm									x
5307-0114	Head, volume 1,4 ml, MicroSampler tube-Ø 12 mm									x
5307-0115	Head, volume 1,5 ml, MicroSampler tube-Ø 12 mm									x
5307-0116	Head, volume 1,6 ml, MicroSampler tube-Ø 12 mm									x
5307-0117	Head, volume 1,7 ml, MicroSampler tube-Ø 12 mm									x
5307-0118	Head, volume 1,8 ml, MicroSampler tube-Ø 12 mm									x
5307-0119	Head, volume 1,9 ml, MicroSampler tube-Ø 12 mm									x
5307-0120	Head, volume 2,0 ml, MicroSampler tube-Ø 12 mm									x
5307-0135	Head, volume 0,35 ml, MicroSampler tube-Ø 12 mm									x
5307-0201	Head, volume 1 ml, MicroSampler tube-Ø 25 mm									x
5307-0202	Head, volume 2 ml, MicroSampler tube-Ø 25 mm									x
5307-0203	Head, volume 3 ml, MicroSampler tube-Ø 25 mm									x
5307-0204	Head, volume 4 ml, MicroSampler tube-Ø 25 mm									x
5307-0205	Head, volume 5 ml, MicroSampler tube-Ø 25 mm									x
5307-0206	Head, volume 6 ml, MicroSampler tube-Ø 25 mm									x
5307-0210	Head, volume 10 ml, MicroSampler tube-Ø 25 mm									x
5307-0211	Head, volume 11 ml, MicroSampler tube-Ø 25 mm									x
5307-0212	Head, volume 12 ml, MicroSampler tube-Ø 25 mm									x
5307-0213	Head, volume 13 ml, MicroSampler tube-Ø 25 mm									x
5307-0214	Head, volume 14 ml, MicroSampler tube-Ø 25 mm									x
5307-0215	Head, volume 15 ml, MicroSampler tube-Ø 25 mm									x
5307-0216	Head, volume 16 ml, MicroSampler tube-Ø 25 mm									x
5307-0217	Head, volume 17 ml, MicroSampler tube-Ø 25 mm									x
5307-0218	Head, volume 18 ml, MicroSampler tube-Ø 25 mm									x
5307-0219	Head, volume 19 ml, MicroSampler tube-Ø 25 mm									x
5307-0220	Head, volume 20 ml, MicroSampler tube-Ø 25 mm									x
5307-0221	Head, volume 0,15 ml, MicroSampler tube-Ø 12 mm									x
5307-0225	Head, volume 2,5 ml, MicroSampler tube-Ø 25 mm									x
5307-0235	Head, volume 3,5 ml, MicroSampler tube-Ø 25 mm									x
5307-1055	MicroSampler, 316L, tube-Ø 12 mm, length 55 cm									x
5307-1085	MicroSampler, 316L, tube-Ø 12 mm, length 85 cm									x
5307-1120	MicroSampler, 316L, tube-Ø 12 mm, length 120 cm									x
5307-2055	MicroSampler, 316L, tube-Ø 25 mm, length 55 cm									x
5307-2085	MicroSampler, 316L, tube-Ø 25 mm, length 85 cm									x
5307-2120	MicroSampler, 316L, tube-Ø 25 mm, length 120 cm									x

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



Article	Description	Alu	PTFE	PA	1.4034	1.4301 304	1.4305	1.4310	1.4571 316Ti	1.4404 316L
5309-0110	Lowering cable AISI304/PTFE, length 10 m, Ø 2 mm		x			x				
5309-0125	Lowering cable AISI304/PTFE, length 25 m, Ø 2 mm		x			x				
5309-0150	Lowering cable AISI304/PTFE, length 50 m, Ø 2 mm		x			x				
5309-0301	Lowering cable 316L, length 10 m, Ø 1,25 mm									x
5309-0302	Lowering cable 316L, length 25 m, Ø 1,25 mm									x
5309-0305	Lowering cable 316L, length 50 m, Ø 1,25 mm									x
5309-0510	Stainless steel chain 316L, length 10 m, Ø 2,3mm									x
5309-0525	Stainless steel chain 316L, length 25 m, Ø 2,3mm									x
5309-0550	Stainless steel chain 316L, length 100 m, Ø 2,3mm									x
5310-0001	Core Sampler, 316L, length 200mm, Ø 15mm		x							x
5310-0002	Core Sampler, 316L, length 400mm, Ø 25mm		x							x
5310-0003	Core Sampler, 316L, length 800mm, Ø 25mm		x							x
5310-0150	All-layer drill, 316L, length 150 cm, volume 260ml									x
5311-0001	QualiRod conical, drill. d. 110mm, EN/ISO 707					x				
5311-0002	QualiRod conical, drill. d. 150mm, EN/ISO 707					x				
5311-0101	QualiRod cylindrical, 304, drilling depth 20 cm					x				
5313-0180	Spatula, blade 1.4034, wooden handle, l. 180 mm					x				
5313-0300	Spatula, blade 1.4034, wooden handle, l. 300 mm					x				
5316-0055	Multi sampler, 316L/PTFE, 55cm, 3 chambers per 14ml		x							x
5316-0085	Multi sampler, 316L/PTFE, 85cm, 3 chambers per 17ml		x							x
5316-0150	Multi sampler, 316L/PTFE, 150cm, 5 chambers/ 17 ml		x							x
5316-1055	Uno sampler, 316L/PTFE, 55cm, 1 chamber 17ml volume		x							x
5316-1085	Uno sampler, 316L/PTFE, 85cm, 1 chamber 17ml volume		x							x
5316-1150	Uno sampler, 316L/PTFE, 150cm, 1 chamber 17 ml vol		x							x
5316-2055	All-layer-sampler, 316L, 55 cm, 3 openings, 70 ml									x
5316-2085	All-layer-sampler, 316L, 85 cm, 3 openings, 130 ml									x
5316-2150	All-layer-sampler, 316L, 150 cm, 5 openings, 260ml									x
5316-3055	Multi sampler, 316L, 55cm, 3 chambers per 14ml vol									x
5316-3085	Multi sampler, 316L, 85cm, 3 chambers per 17ml vol									x
5316-3150	Multi sampler, 316L, 150cm, 5 chambers per 17 ml									x
5317-0015	Depth stop, inner dia. 25 mm, AISI 316L									x
5317-0020	Depth stop, inner dia. 50 mm, aluminium	x				x				x
5317-0030	Emptying device, AISI304, length 100 cm, complete					x				
5317-0033	Emptying funnel separate, 304					x				
5317-1058	Multi sampler, Alu, 55 cm, 3 chambers per 14 ml	x								
5317-1085	Multi sampler, Alu, 85 cm, 3 chambers per 17 ml	x								
5317-1150	Multi sampler, Alu, 150cm, 5 chambers per 17 ml	x								

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



Article	Description	Alu	PTFE	PA	1.4034	1.4301 304	1.4305	1.4310	1.4571 316Ti	1.4404 316L
5317-2055	All-layer-sampler, Alu, 55 cm, 3 openings, 70 ml	x								
5317-2085	All-layer-sampler, Alu, 85 cm, 3 openings, 130 ml	x								
5317-2150	All-layer-sampler, Alu, 150 cm, 5 openings, 260ml	x								
5318-0145	Zone sampler Spiralus, Alu/PA, 145cm, 770ml volume	x		x						
5318-0200	Zone sampler Spiralus, Alu/PA, 200cm, 1100 ml vol.	x		x						
5319-1085	Jumbo sampler, all layer, Alu, 85cm, 880 ml volume	x								
5319-1150	Jumbo sampler, all layer, Alu, 150 cm, 1700ml vol.	x								
5319-1250	Jumbo sampler, all layer, Alu, 250 cm, 2900ml vol.	x								
5319-2085	Jumbo sampler, closed chambers, Alu/PTFE, 85 cm	x	x							
5319-2150	Jumbo sampler, closed chambers, Alu/PTFE, 150 cm	x	x							
5319-2250	Jumbo sampler, closed chambers, Alu/PTFE, 250 cm	x	x							
5320-0150	Sample spoon, AISI 304, 150 mm, 2 ml, autoclavable					x				
5320-0180	Sample spoon, AISI 304, 180 mm, 2 ml, autoclavable					x				
5320-0210	Sample spoon, AISI 304, 210 mm, 2 ml, autoclavable					x				
5320-0250	Sample spoon, AISI 304, 250 mm, 9 ml, autoclavable					x				
5320-0300	Sample-spoon, V2A, 300 mm, 6 ml, autoclavable					x				
5320-0400	Sample spoon, AISI 304, 400 mm, 9 ml, autoclavable					x				
5320-0500	Sample spoon, AISI 304, 500 mm, 9 ml, autoclavable					x				
5321-0300	Sampling lance Milky, 316L, 480 mm, 150 ml, Ø 28mm									x
5322-1300	QuickPicker 316L, 300 mm, 75 ml volume, ISTA-Norm									x
5323-0190	Ice borer, 316L, inner-xouter-Ø 8x10mm, depth 10cm								x	x
5323-0200	Ice borer, 316L, inner-xouter-Ø 16x21mm,depth 20cm								x	x
5323-2010	Ice borer, 316L, inner-xouter-Ø 50x54mm,depth 20cm								x	x
5324-0001	PharmaSpoon, 316L, length 650 mm, hand polished									x
5324-0010	PharmaScoop, 316L, 10 ml vol., Ø 25 mm,length 45mm									x
5324-0011	PharmaScoop, 316L, 50 ml vol., Ø 38 mm,length 70mm									x
5324-0012	PharmaScoop, 316L, 100ml vol., Ø 50mm,length 100mm									x
5324-0013	PharmaScoop, 316L, 750ml vol., Ø 95mm,length 210mm									x
5324-0014	PharmaScoop, 316L, 1000 ml vol., Ø110 mm, 235 mm									x
5324-0015	PharmaScoop, 316L, 1250 ml vol., Ø125 mm, 195 mm									x
5324-0016	PharmaScoop, 316L, 1700 ml vol., Ø135 mm, 230 mm									x
5324-0017	PharmaScoop, 316L, 3000ml vol., Ø180mm, l. 215mm									x
5324-0018	PharmaScoop, 316L, 500ml vol., Ø 95mm,length 140mm									x
5324-0019	PharmaScoop with long handle, 316L, 1000 ml									x
5326-0001	Liquid-Sampler, open with thumb, 316L, 50 ml									x
5326-0002	Liquid-Sampler, open with thumb, 316L, 100 ml									x
5326-0003	Liquid-Sampler, place onto bottle, 316L, 50 ml									x

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



Article	Description	Alu	PTFE	PA	1.4034	1.4301 304	1.4305	1.4310	1.4571 316Ti	1.4404 316L
5326-0004	Liquid-Sampler, place onto bottle, 316L, 100 ml									x
5334-8000	Stainless steel measuring spoons set (1,25-15 ml)					x				
5337-2020	Sampling drill for stalk material, Ø25mm, 200mm				x	x				x
5337-2040	Sampling drill for stalk material, Ø25mm, 400mm				x	x				x
5337-4020	Sampling drill for stalk material, Ø45mm, 200mm	x			x	x				
5337-4040	Sampling drill for stalk material, Ø45mm, 400mm	x			x	x				
5338-0410	StickProof, 316L, length 410 mm, Ø 25mm, ISTA-Norm									x
5338-1100	StickProof, 316L, length 1100mm, Ø 25mm, ISTA-Norm									x
5339-0300	IceSampler, 316L, length 300mm, vol. 50ml, Ø 38mm									x
5339-0550	IceSampler, 316L, length 550mm, vol. 50ml, Ø 38mm									x
5339-1050	IceSampler, 316L, length 1050mm, vol. 50ml, Ø 38mm									x
5340-1100	SiloPicker, 316L, vol. 100 ml, 130 cm, outer-Ø50mm									x
5341-1002	Suction head, suction drill 2mm, Mini ViscoSampler									x
5341-1004	Suction head, suction drill 4mm, Mini ViscoSampler									x
5341-1006	Suction head, suction drill 6mm, Mini ViscoSampler									x
5341-1060	Mini ViscoSampler, 316L, 60 cm, vol. 30 ml, Ø 15mm		x							x
5341-1100	Mini ViscoSampler, 316L, 100 cm, vol. 50 ml, Ø15mm		x							x
5341-1150	Mini ViscoSampler, 316L, 150 cm, vol. 75 ml, Ø15mm		x							x
5346-0010	PharmaPicker case, complete in transport box									x
5346-1000	PharmaPicker, 316L, w/collection cylinder 0,1-1,2ml									x
5346-1001	Collection cylinder for 0,1-1,2ml, w/o head									x
5346-1010	Volume head, capacity 0,1 ml, PharmaPicker									x
5346-1020	Volume head, capacity 0,2 ml, PharmaPicker									x
5346-1060	Volume head, capacity 0,6 ml, PharmaPicker									x
5346-1100	Volume head, capacity 1 ml, PharmaPicker									x
5346-2000	PharmaPicker, 316L, w/coll.cylinder for 1,25-2,5ml									x
5346-2001	Collection cylinder for 1,25-2,5ml, w/o head									x
5346-2200	Volume head, capacity 2 ml, PharmaPicker									x
5346-2250	Volume head, capacity 2,5 ml, PharmaPicker									x
5346-3010	Powder sampler Mini, AISI316L, Ø6mm, 300mm, 0,01ml								x	x
5346-3050	Powder sampler Mini, AISI316L, Ø6mm, 300mm, 0,05ml								x	x
5346-3100	Powder sampler Mini, AISI316L, Ø6mm, 300mm, 0,1ml								x	x
5347-0005	Scratcher, 316L, 5 ml, flute length 22 mm, Ø 25 mm									x
5347-0010	Scratcher, 316L, 10 ml, flute length 45 mm, Ø 25mm									x
5347-0025	Scratcher, 316L, 25 ml, flute length 110 mm, Ø25mm									x
5351-1000	All-layer-scratcher, 316L, Ø 25mm, 100cm, vol.180ml		x			x				x
5354-1000	Stainless steel scoop, 304, 115cm, capacity 1000ml					x				

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



Article	Description	Alu	PTFE	PA	1.4034	1.4301 304	1.4305	1.4310	1.4571 316Ti	1.4404 316L
5354-2000	Stainless steel beaker, 304, 1000 ml, TeleScoop					x				
5354-2010	Stainless steel beaker, extended, 1000ml, TeleScoop					x				
5354-3000	Pendulum beaker 304, capacity 1000 ml, TeleScoop					x	x			
5356-1500	EasyScoop, aluminium, 1500 ml, inner grip	x								
5356-2000	EasyScoop, aluminium, 2000 ml, inner grip	x								
5357-0001	Sampler Liquid CupSampler, 316L, 100 ml, 100 cm									x
5358-0000	Sample insert, 0,0 ml, Novartos Multi									x
5358-0005	Sample insert, 0,5 ml, Novartos Multi									x
5358-0010	Sample insert, 1 ml, Novartos Multi									x
5358-0015	Sample insert, 1,5 ml, Novartos Multi									x
5358-0020	Sample insert, 2 ml, Novartos Multi									x
5358-0100	Sample insert, 0,0 ml, Novartos Uno									x
5358-0105	Sample insert, 0,5 ml, Novartos Uno									x
5358-0110	Sample insert, 1 ml, Novartos Uno									x
5358-0120	Sample insert, 2 ml, Novartos Uno									x
5358-0130	Sample insert, 3 ml, Novartos Uno									x
5358-0140	Sample insert, 4 ml, Novartos Uno									x
5358-0600	Novartos Uno, 60cm, double, 2 chambers, w/o inserts									x
5358-0602	Novartos Multi, 60cm, double, 8 chambers, w/o inserts									x
5358-0603	Novartos Multi, 60cm, triple, 12 chambers, w/o inserts									x
5358-0630	Novartos Uno, 60cm, triple, 3 chambers, w/o inserts									x
5358-1000	Novartos Multi, 100cm, double, 14 chambers, w/o inser.									x
5358-1001	Novartos Uno, 100cm, double, 2 chambers, w/o inser.									x
5358-1003	Novartos Multi, 100cm, triple, 21 chambers, w/o inser.									x
5358-1011	Novartos 316L, head									x
5358-1030	Novartos Uno, 100cm, triple, 3 chambers, w/o inser.									x
5358-1400	Novartos Uno, 140cm, double, 2 chambers, w/o inser.									x
5358-1402	Novartos Multi, 140cm, double, 18 chambers, w/o inser.									x
5358-1403	Novartos Multi, 140cm, triple, 27 chambers, w/o inser.									x
5358-1430	Novartos Uno, 140cm, triple, 3 chambers, w/o inser.									x
5359-1600	Valve lance, 316L, length 550 mm, Ø 32 mm									x
5359-1800	Valve lance, 316L, length 800 mm, Ø 32 mm									x
5360-0600	EasySampler, 316L, length 100 cm, volume 60 ml									x
5360-1000	EasySampler, 316L, length 104,5 cm, volume 200 ml									x
5369-0300	Pigment lance, flute open. 25%, length 300mm, 145ml									x
5369-0310	Pigment lance, flute open. 50%, flute 300mm, 85 ml									x
5369-0600	Pigment lance, flute open. 25%, flute 600mm, 312ml									x

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



Article	Description	Alu	PTFE	PA	1.4034	1.4301 304	1.4305	1.4310	1.4571 316Ti	1.4404 316L
5369-0610	Pigment lance, flute open. 50%, flute 600mm, 179ml									x
5369-0900	Pigment lance, flute open. 25%, flute 900mm, 480ml									x
5369-0910	Pigment lance, flute open. 50%, flute 900mm, 273ml									x
5370-0120	Scoop AISI AISI304, capacity 120 ml, Ø 50 mm					x				
5370-0450	Scoop AISI AISI304, capacity 450 ml, Ø 75 mm					x				
5370-1000	Scoop AISI AISI304, capacity 1000 ml, Ø 100 mm					x				
5370-1500	Scoop AISI AISI304, capacity 1500 ml, Ø 120 mm					x				
5370-2000	Scoop AISI AISI304, capacity 2500 ml, Ø 150 mm					x				
5371-0100	FoodScoop, 304, 100ml, total 215mm, LxW 120x60mm					x				
5371-0200	FoodScoop, 304, 200ml, total 230mm, LxW 140x75mm					x				
5371-0250	FoodScoop, 304, 250ml, total 245mm, LxW 155x85mm					x				
5371-0350	FoodScoop, 304, 350ml, total 290mm, LxW 180x100mm					x				
5371-0500	FoodScoop, 304, 500ml, total 350mm, LxW 220x130mm					x				
5371-1000	FoodScoop, 304, 1000ml, total 395mm, LxW 260x142mm					x				
5371-1500	FoodScoop, 304, 1500ml, total 425mm, LxW 290x160mm					x				
5371-2000	FoodScoop, 304, 2000ml, total 465mm, LxW 328x177mm					x				
5371-5002	Ladle, stainless steel AISI 304, 200 ml					x				
5371-5004	Ladle, stainless steel AISI 304, 620 ml					x				
5380-0506	Extension rod, 500 mm, M6 thread, PharmaPicker									x
5380-0508	Extension rod, 500 mm, M8 thread, SiloPicker									x
5380-1006	Extension rod, 1000 mm, M6 thread, PharmaPicker									x
5380-1008	Extension rod, 1000 mm, M8 thread, SiloPicker									x
5382-0040	Sampler Tubus, 316L, 40 cm, Ø 40 mm, ISTA standard									x
5382-0080	Sampler Tubus, 316L, 80 cm, Ø 40 mm, ISTA standard									x
5383-1000	Spatula stainless steel, 316L, 300 mm									x
5386-0001	Micro spatula stainless steel, 304, length 170 mm					x				
5386-0010	Palette knife spatula stainl. steel, LxW 200x21 mm					x				
5386-0012	Palette knife spatula stainl. steel, LxW 250x25 mm					x				
5386-0020	Spoon spatula stainless steel, AISI304, 120mm, 2ml					x				
5386-0022	Spoon spatula stainless steel, AISI304, 180mm, 4ml					x				
5386-0024	Spoon spatula stainless steel, AISI 304, 250mm, 7ml					x				
5386-0025	Spoon spatula stainless steel, AISI304, 300mm, 11ml					x				
5386-0030	Scraper stainless steel, AISI 316, LxW 230 x 75mm									x
5386-0035	Scraper stainless steel, AISI 316, LxW 270 x 100mm									x
5386-0041	Double spatula square/rounded, AISI 304, 150 mm					x				
5386-0043	Double spatula square/rounded, AISI 304, 210 mm					x				
5386-0045	Double spatula square/rounded, AISI 304, 300 mm					x				

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



Article	Description	Alu	PTFE	PA	1.4034	1.4301 304	1.4305	1.4310	1.4571 316Ti	1.4404 316L
5386-0048	Double spatula square/rounded, AISI 304, 500 mm					x				
5386-0052	Micro double spatula square/rounded, 304, 150x2 mm					x				
5386-0055	Micro double spatula square/rounded, 304, 150x3 mm					x				
5386-0058	Micro double spatula square/rounded, 304, 150x4 mm					x				
5386-0064	Micro double spatula square/rounded, 304, 210x5 mm					x				
5386-0067	Micro double spatula square/rounded, 304, 210x6 mm					x				
5386-0070	Micro double spatula square/rounded, 304, 210x9 mm					x				
5386-0080	Chattaway double spatula square/round bent, 125 mm					x				
5386-0081	Chattaway double spatula square/round bent, 150 mm					x				
5386-0082	Chattaway double spatula square/round bent, 180 mm					x				
5386-0083	Chattaway double spatula square/round bent, 200 mm					x				
5386-0084	Chattaway double spatula square/round bent, 235 mm					x				
5386-0091	Powder spatula square/half-shell, AISI 304, 150 mm					x				
5386-0093	Powder spatula square/half-shell, AISI 304, 210 mm					x				
5386-0096	Micro powder spatula square/half-shell, 150 mm					x				
5386-0098	Micro powder spatula square/half-shell, 210 mm					x				
5386-0100	Forceps, stainless steel, sharp, bent form, 105 mm					x				
5386-0102	Forceps, stainless steel, sharp, bent form, 130 mm					x				
5386-0104	Forceps, stainless steel, sharp, bent form, 160 mm					x				
5386-0200	Forceps, stainless steel, sharp, straight, 105 mm					x				
5386-0202	Forceps, stainless steel, sharp, straight, 130 mm					x				
5386-0204	Forceps, stainless steel, sharp, straight, 160 mm					x				
5386-0300	Forceps, stainless steel, blunt, straight, 105 mm					x				
5386-0302	Forceps, stainless steel, blunt, straight, 130 mm					x				
5386-0304	Forceps, stainless steel, blunt, straight, 160 mm					x				
5386-0401	Cover glass tweezer, AISI 304, 105 mm					x				
5386-0410	Cover glass tweezer, curved, AISI 304, 105 mm					x				
5386-1001	Double spatula, flat, square, AISI 304, 100 mm					x				
5386-1002	Double spatula, flat, square, AISI 304, 150 mm					x				
5386-1012	Micro spoon spatula, square, 304, 150 mm, 0.1 ml					x				
5386-1012	Micro spoon spatula, square, 304, 150 mm, 0.1 ml					x				
5386-1013	Micro spoon spatula, square, 304, 150 mm, 0.07 ml					x				
5386-1013	Micro spoon spatula, square, 304, 150 mm, 0.07 ml					x				
5386-1014	Micro spoon spatula, square, 304, 150 mm, 0.04 ml					x				
5386-1014	Micro spoon spatula, square, 304, 150 mm, 0.04 ml					x				
5386-1015	Micro spoon spatula, square, 304, 200 mm, 0.12 ml					x				
5386-1015	Micro spoon spatula, square, 304, 200 mm, 0.12 ml					x				

ZERTIFIKAT CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



Article	Description	Alu	PTFE	PA	1.4034	1.4301 304	1.4305	1.4310	1.4571 316Ti	1.4404 316L
5386-1021	Spoon spatula, round handle, AISI304, 150 mm, 2 ml					x				
5386-1023	Spoon spatula, round handle, 304, 210 mm, 3.5 ml					x				
5386-1025	Spoon spatula, round handle, AISI304, 300 mm, 5 ml					x				
5386-1031	Chemical spoon, round handle, AISI304, 150 mm, 2ml					x				
5386-1033	Chemical spoon, round handle, 304, 210 mm, 3.5 ml					x				
5386-1035	Chemical spoon, round handle, AISI304, 300 mm, 5ml					x				
5386-1041	Laboratory spoon, double-ended, 304, 150mm,2/1.5ml					x				
5386-1043	Laboratory spoon, double-ended, 304, 210mm,3.5/2ml					x				
5386-1045	Laboratory spoon, double-ended, 304, 300mm,5/3.5ml					x				
5386-4000	Beaker tongs, AISI 304, silicone coating, 250 mm					x				
5388-1000	BeefSteaker, 304, Ø 25mm				x					
5388-1010	Drill head for carcass sampling				x					
6213-0085	Aluminium container with lid, 85 ml	x								
6213-0175	Aluminium container with lid, 175 ml	x								
6213-0375	Aluminium container with lid, 375 ml	x								
7108-0003	Measuring cup AISI 304 with handle/spout, 250 ml					x				
7108-0005	Measuring cup AISI 304 with handle/spout, 500 ml					x				
7108-0010	Measuring cup AISI 304 with handle/spout, 1000 ml					x				
7108-0020	Measuring cup AISI 304 with handle/spout, 2000 ml					x				
7121-0100	Laboratory beaker stainless steel, AISI 304, 100ml					x				
7121-0250	Laboratory beaker stainless steel, AISI 304, 250ml					x				
7121-0500	Laboratory beaker stainless steel, AISI 304, 500ml					x				
7121-1000	Laboratory beaker stainless steel, AISI 304,1000ml					x				
7121-3000	Laboratory beaker stainless steel, AISI 304,3000ml					x				
7121-5000	Laboratory beaker stainless steel, AISI 304,5000ml					x				
9604-0001	Funnel, stainless steel, AISI 304, outer-Ø 60 mm					x				
9604-0002	Funnel, stainless steel, AISI 304, outer-Ø 80 mm					x				
9604-0003	Funnel, stainless steel, AISI 304, outer-Ø 100 mm					x				
9604-0004	Funnel, stainless steel, AISI 304, outer-Ø 120 mm					x				
9604-0005	Funnel, stainless steel, AISI 304, outer-Ø 150 mm					x				
9604-0006	Funnel, stainless steel, AISI 304, outer-Ø 200 mm					x				
9604-0007	Funnel, stainless steel, AISI 304, outer-Ø 250 mm					x				
9604-0008	Funnel, stainless steel, AISI 304, outer-Ø 300 mm					x				
9604-0100	Powder funnel, AISI 304, outer Ø 90 mm					x				
9635-1001	Test tube stand, AISI 304, 13x12 mm					x				
9635-1002	Test tube stand, AISI 304, 13x24 mm					x				
9635-1005	Test tube stand, AISI 304, 16x12 mm					x				

ZERTIFIKAT CERTIFICATE

Pumpen | Probennehmer | Laborbedarf
für Labor, Industrie und Wissenschaft
Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science



Article	Description	Alu	PTFE	PA	1.4034	1.4301 304	1.4305	1.4310	1.4571 316Ti	1.4404 316L
9635-1006	Test tube stand, AISI 304, 16x24 mm					x				
9635-1009	Test tube stand, AISI 304, 20x12 mm					x				
9635-1010	Test tube stand, AISI 304, 20x24 mm					x				
9635-2001	Test tube holder Z-shaped, AISI 304, 85x85 mm					x				
9635-2002	Test tube holder Z-shaped, AISI 304, 85x165 mm					x				
9635-2003	Test tube holder Z-shaped, AISI 304, 110x110 mm					x				
9635-2004	Test tube holder Z-shaped, AISI 304, 110x220 mm					x				