

5202-x000 PumpMaster
5202-x001 Discharge hose w/ stopcock, PumpMaster

Herewith we confirm that parts coming in contact with the medium are made of the on page 2 listed materials. The pumps are not suitable for hydrochloric acid (stainless steel spring).

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of certificate is 2284-08.

14. April 2021



Bürkle GmbH, Bad Bellingen,
Martin Saint-Denis, Managing Director



Article	Description	PP	NBR	EPDM	FKM
5202-1000	PumpMaster acids/ chemical liquids, FKM, green				
	pump body	x			
	tube	x			
	gasket				x
5202-2000	PumpMaster non-aggressive liquids, PP/EPDM, red				
	pump body	x			
	tube	x			
	gasket			x	
5202-3000	PumpMaster petrochem. liquids, PP/Nitrile r., blue				
	pump body	x			
	tube	x			
	gasket		x		

Article	Description	PP	NBR	EPDM	FKM
5202-1001	Discharge hose w/ stopcock, PumpMaster acids				
	hose				x
	stopcock	x			
	gasket				x
5202-2001	Discharge hose w/ stopcock, PumpMaster liquids				
	hose			x	
	stopcock	x			
	gasket			x	
5202-3001	Discharge hose w/ stopcock, PumpMaster petrochem.				
	hose		x		
	stopcock	x			
	gasket		x		