

Declaration of Compliance

Our supplier confirms that the below mentioned products comply with:

- (EC) No. 1935/2004
- (EU) No. 10/2011
- (EU) 2018/213
- (EC) No. 2023/2006
- 94/62/EC
- free of Bisphenol A

The results of global and specific migration are under the legal limit. The user shall verify himself that the product is suitable for the intended food to be packed beyond the stipulations of the Directives. In particular it is emphasized that in case of printed material, no contact is allowed between the printing ink and food product.

Further our supplier confirms that the below mentioned products comply with:

- EP 3.1.3 (except for 035x-, 0304-xxxx)
- free of Phthalates (except for 035x-xxxx)
- free of BSE (except for 035x-xxxx)
- FDA 21 §177.1520 (except for 035x-xxxx)
- 2011/65/EU ROHS (except for 0357-xxxx)

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of certificate is 2284-08

21. March 2019



Bürkle GmbH, Bad Bellingen,
Martin Saint-Denis, Managing Director



Item no.	Description
0301-0100	Multi-Bottle w/ graduation, HDPE, 100 ml, w/ cap
0301-0250	Multi-Bottle w/ graduation, HDPE, 250 ml, w/ cap
0301-0500	Multi-Bottle w/ graduation, HDPE, 500 ml, w/ cap
0301-1000	Multi-Bottle w/ graduation, HDPE, 1000 ml, w/ cap
0302-0200	Narrow-necked bottle, LDPE transp., 200 ml, w/ cap
0302-3000	Narrow-necked bottle, LDPE transp., 3000 ml, w/cap
0302-5000	Narrow-necked bottle, LDPE transp., 5000 ml, w/cap
0304-0500	UN-bottle, HDPE, 500ml, w/ tamper-evident closure
0304-1000	UN-bottle, HDPE, 1000ml, w/ tamper-evident closure
0304-1002	UN-Flasche, HDPE, 1000ml, o. Verschluss schwarz
0304-2500	UN-bottle, HDPE, 2500ml, w/ tamper-evident closure
0308-0050	Round bottle 50ml
0308-0100	Round bottle w/ snap closure, HDPE, 100 ml
0308-0250	Round bottle w/ snap closure, HDPE, 250 ml
0308-0500	Round bottle w/ snap closure, HDPE, 500 ml
0308-1000	Round bottle w/ snap closure, HDPE, 1000 ml
0317-2300	Wide-mouth container w/handle, HDPE, 2300ml, w/cap
0317-4400	Wide-mouth container w/handle, HDPE, 4400ml, w/cap
0318-0200	Wide-necked bottle, LPDE, round, 200 ml, w/ cap
0318-0300	Wide mouth bottle 300 ml
0318-0750	Wide-necked bottle, LPDE, round, 750 ml, w/ cap
0318-1500	Wide-necked bottle, LPDE, round, 1500 ml, w/ cap
0338-0250	Wide-necked square bottle, HDPE, 250 ml, w/ cap
0338-0500	Wide-necked square bottle, HDPE, 500 ml, w/ cap
0340-0100	Narrow-necked reagent bottle, HDPE, 100 ml, w/ cap
0340-0250	Narrow-necked reagent bottle, HDPE, 250 ml, w/ cap
0340-0500	Narrow-necked reagent bottle, HDPE, 500 ml, w/ cap
0340-1000	Narrow-necked reagent bottle, HDPE, 1000ml, w/ cap
0340-2500	Narrow-necked reagent bottle, HDPE, 2500ml, w/ cap
0342-0050	Wide-necked reagent bottle, HDPE, 50 ml, w/ cap
0342-0100	Wide-necked reagent bottle, HDPE, 100 ml, w/ cap
0342-0250	Wide-necked reagent bottle, HDPE, 250 ml, w/ cap
0342-0500	Wide-necked reagent bottle, HDPE, 500 ml, w/ cap
0342-0750	Wide-necked reagent bottle, HDPE, 750 ml, w/ cap
0342-1000	Wide-necked reagent bottle, HDPE, 1000 ml, w/ cap
0342-1500	Wide-necked reagent bottle, HDPE, 1500 ml, w/ cap
0342-2500	Wide-necked reagent bottle, HDPE, 2500 ml, w/ cap
0342-4000	Wide-necked reagent bottle, HDPE, 4000 ml, w/ cap
0355-0050	Wide-necked box, square, PVC transp., 50ml, w/ cap
0355-0100	Wide-necked box, square, PVC transp., 100ml, w/cap
0355-0111	Weithalsflasche, PVC, 100ml, ohne Verschluss
0355-0200	Wide-necked box, square, PVC transp., 200ml, w/cap
0355-0211	Weithalsdose, PVC, 200ml, ohne Verschluss
0355-0300	Wide-necked box, square, PVC transp., 300ml, w/cap
0355-0500	Wide-necked box, square, PVC transp., 500ml, w/cap
0355-0750	Wide-necked box, square, PVC transp., 750ml, w/cap
0355-1000	Wide-necked box square, PVC transp., 1000ml, w/cap
0355-2000	Wide-necked box square, PVC transp., 2000ml, w/cap

Item no.	Description
0356-0100	Wide-necked box, square, PVC amber, 100 ml, w/ cap
0356-0200	Wide-necked box, square, PVC amber, 200 ml, w/ cap
0356-0100	Wide-necked box, square, PVC amber, 100 ml, w/ cap
0356-0200	Wide-necked box, square, PVC amber, 200 ml, w/ cap
0356-1000	Wide neck box amber 1000 ml
0357-0050	Wide-necked box, square, PETG clear, 50 ml, w/ cap
0357-0100	Wide-necked box, square, PETG clear, 100 ml, w/cap
0357-0200	Wide-necked box, square, PETG clear, 200 ml, w/cap
0357-0500	Wide-necked box, square, PETG clear, 500 ml, w/cap
0357-1000	Wide-necked box, square, PETG clear, 1000ml, w/cap
0357-2000	Wide neck rect. box 2000 ml
0401-0005	Storage bottle w/o thread. con., HDPE, 5 l, w/ cap
0401-0010	Storage bottle w/o thread. con., HDPE, 10l, w/ cap
0401-0025	Storage bottle w/o thread. con., HDPE, 25l, w/ cap
0401-0060	Storage bottle w/o thread. con., HDPE, 60l, w/ cap
0402-0005	Storage bottle w/ thread. con., HDPE, 5 l, w/ cap
0402-0010	Storage bottle w/ thread. con., HDPE, 10 l, w/ cap
0402-0025	Storage bottle w/ thread. con., HDPE, 25 l, w/ cap
0402-0060	Storage bottle w/ thread. con., HDPE, 60 l, w/ cap
1406-1005	Wide-mouth container, HDPE, 5 l, w/ cap
1406-1010	Wide-mouth container, HDPE, 10 l, w/ cap
2306-0010	Laboratory bucket, PE white, w/ spout, 10 l
2306-0015	Laboratory bucket, PE white, w/ spout, 15 l
6282-0100	Wide-necked box round, HDPE, 100ml, Ø 48 mm, w/cap
6282-0250	Wide-necked box round, HDPE, 250ml, Ø 65 mm, w/cap
6282-0500	Wide-necked box round, HDPE, 500ml, Ø 80 mm, w/cap
6282-0750	Wide-necked box round, HDPE, 750ml, Ø 94 mm, w/cap
6282-1000	Wide-necked box round, HDPE, 1000ml, Ø100mm, w/cap
6282-1500	Wide-necked box round, HDPE, 1500ml, Ø114mm, w/cap
6282-2000	Wide-necked box round, HDPE, 2000ml, Ø126mm, w/cap