

Material certificate

The articles on the following pages are manufactured from stainless steel (V2A or V4A) or Aluminum, some in combination with PTFE. They are all designed to have a smooth finish and to be internally crevice free. They can be easily cleaned in line with GMP.

The declaration for stainless steel is based on REGULATION (EC) No 1935/2004 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL "The principle underlying this regulation is that any material or article intended to come into contact, directly or indirectly, with food must be sufficiently inert to preclude substances from being transferred to food in quantities large enough to endanger human health or to bring about an unacceptable change in the composition of the food or a deterioration in its organoleptic properties."

Stainless steel is in general widely used and suitable for most foodstuffs. If using highly corrosive or aggressive materials, please check our List of "Chemical-Resistance" :

https://www.buerkle.de/files_pdf/wissenswertes/chemical_resistance_en.pdf or contact us for suitability.

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up.

The supplier - Bürkle GmbH in Bad Bellingen/Germany - is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of the certificate is 2284-08.

24. January 2019



Bürkle GmbH, Bad Bellingen,
Martin Saint-Denis, Managing Director



Article	Description	V2A	V4A	Alu	PTFE	PA
2309-0010	Bucket, stainless steel, V2A, capacity 10 l	X				
2309-0015	Bucket, stainless steel, V2A, capacity 15 l	X				
2601-1001	Safety can V4A, dosage spout, relief valve, 1 l		X		X	
2601-1002	Safety can V4A, dosage spout, relief valve, 2 l		X		X	
2601-1005	Safety can V4A, dosage spout, relief valve, 5 l		X		X	
2601-2010	Safety can V4A, dosage spout, relief valve, 10 l		X		X	
2601-2020	Safety can V2A, dosage spout, relief valve, 20 l	X			X	
2601-3005	Safety jerrycan V4A, screw cap, relief valve, 5 l		X		X	
2601-3010	Safety jerrycan V4A, screw cap, relief valve, 10 l		X		X	
2601-3020	Safety jerrycan V2A, screw cap, relief valve, 20 l	X			X	
2601-4010	Safety storage cont., V4A, tab, relief valve, 10 l		X		X	
2601-4025	Safety storage cont., V2A, tab, relief valve, 25 l	X			X	
2601-5050	Safety storage container V2A, relief valve, 50 l	X			X	
2601-5125	Safety storage container, V2A, relief valve, 25 l	X			X	
2601-6010	Safety storage cont., V4A, tab, level indic., 10 l		X		X	
2601-7010	Safety jerrycan V4A, w/o relief valve, 10 l		X		X	
5301-0500	Extension V2A, 100 cm, SiloDrill	X				
5301-1000	SiloDrill, alu, 150 cm, spiral screw Ø 90mm, 400ml			X		
5301-2000	SiloDrill, V2A, 150 cm, spiral screw Ø 90mm, 400ml	X				
5307-0101	Head, volume 0,1 ml, MicroSampler tube-Ø 12 mm		X			
5307-0102	Head, volume 0,2 ml, MicroSampler tube-Ø 12 mm		X			
5307-0103	Head, volume 0,3 ml, MicroSampler tube-Ø 12 mm		X			
5307-0104	Head, volume 0,4 ml, MicroSampler tube-Ø 12 mm		X			
5307-0105	Head, volume 0,5 ml, MicroSampler tube-Ø 12 mm		X			
5307-0106	Head, volume 0,6 ml, MicroSampler tube-Ø 12 mm		X			
5307-0107	Head, volume 0,7 ml, MicroSampler tube-Ø 12 mm		X			
5307-0108	Head, volume 0,8 ml, MicroSampler tube-Ø 12 mm		X			
5307-0109	Head, volume 0,9 ml, MicroSampler tube-Ø 12 mm		X			
5307-0110	Head, volume 1 ml, MicroSampler tube-Ø 12 mm		X			
5307-0111	Head, volume 1,1 ml, MicroSampler tube-Ø 12 mm		X			
5307-0112	Head, volume 1,2 ml, MicroSampler tube-Ø 12 mm		X			
5307-0113	Head, volume 1,3 ml, MicroSampler tube-Ø 12 mm		X			
5307-0114	Head, volume 1,4 ml, MicroSampler tube-Ø 12 mm		X			
5307-0115	Head, volume 1,5 ml, MicroSampler tube-Ø 12 mm		X			
5307-0116	Head, volume 1,6 ml, MicroSampler tube-Ø 12 mm		X			
5307-0117	Head, volume 1,7 ml, MicroSampler tube-Ø 12 mm		X			
5307-0118	Head, volume 1,8 ml, MicroSampler tube-Ø 12 mm		X			
5307-0119	Head, volume 1,9 ml, MicroSampler tube-Ø 12 mm		X			
5307-0120	Head, volume 2,0 ml, MicroSampler tube-Ø 12 mm		X			
5307-0135	Head, volume 0,35 ml, MicroSampler tube-Ø 12 mm		X			
5307-0201	Head, volume 1 ml, MicroSampler tube-Ø 25 mm		X			
5307-0202	Head, volume 2 ml, MicroSampler tube-Ø 25 mm		X			
5307-0203	Head, volume 3 ml, MicroSampler tube-Ø 25 mm		X			
5307-0204	Head, volume 4 ml, MicroSampler tube-Ø 25 mm		X			
5307-0205	Head, volume 5 ml, MicroSampler tube-Ø 25 mm		X			
5307-0206	Head, volume 6 ml, MicroSampler tube-Ø 25 mm		X			
5307-0210	Head, volume 10 ml, MicroSampler tube-Ø 25 mm		X			
5307-0211	Head, volume 11 ml, MicroSampler tube-Ø 25 mm		X			
5307-0212	Head, volume 12 ml, MicroSampler tube-Ø 25 mm		X			

Article	Description	V2A	V4A	Alu	PTFE	PA
5307-0213	Head, volume 13 ml, MicroSampler tube-Ø 25 mm		X			
5307-0214	Head, volume 14 ml, MicroSampler tube-Ø 25 mm		X			
5307-0215	Head, volume 15 ml, MicroSampler tube-Ø 25 mm		X			
5307-0216	Head, volume 16 ml, MicroSampler tube-Ø 25 mm		X			
5307-0217	Head, volume 17 ml, MicroSampler tube-Ø 25 mm		X			
5307-0218	Head, volume 18 ml, MicroSampler tube-Ø 25 mm		X			
5307-0219	Head, volume 19 ml, MicroSampler tube-Ø 25 mm		X			
5307-0220	Head, volume 20 ml, MicroSampler tube-Ø 25 mm		X			
5307-0221	Head, volume 0,15 ml, MicroSampler tube-Ø 12 mm		X			
5307-0225	Head, volume 2,5 ml, MicroSampler tube-Ø 25 mm		X			
5307-0235	Head, volume 3,5 ml, MicroSampler tube-Ø 25 mm		X			
5307-1055	MicroSampler, V4A, tube-Ø 12 mm, length 55 cm		X			
5307-1085	MicroSampler, V4A, tube-Ø 12 mm, length 85 cm		X			
5307-1120	MicroSampler, V4A, tube-Ø 12 mm, length 120 cm		X			
5307-2055	MicroSampler, V4A, tube-Ø 25 mm, length 55 cm		X			
5307-2085	MicroSampler, V4A, tube-Ø 25 mm, length 85 cm		X			
5307-2120	MicroSampler, V4A, tube-Ø 25 mm, length 120 cm		X			
5310-0150	All-layer drill, V4A, length 150 cm, volume 260 ml		X			
5311-0001	QualiRod conical, V2A, drill. d. 120mm, EN/ISO 707	X				
5311-0002	QualiRod conical, V2A, drill. d. 170mm, EN/ISO 707	X				
5311-0101	QualiRod cylindrical, V2A, drilling depth 20 cm	X				
5313-0180	Spatula, blade V2A, wooden handle, length 180 mm	X				
5313-0300	Spatula, blade V2A, wooden handle, length 300 mm	X				
5316-0055	Multi sampler, V4A/PTFE, 55cm, 3 chambers per 14ml		X		X	
5316-0085	Multi sampler, V4A/PTFE, 85cm, 3 chambers per 17ml		X		X	
5316-0150	Multi sampler, V4A/PTFE, 150cm, 5 chambers/ 17 ml		X		X	
5316-1055	Uno sampler, V4A/PTFE, 55cm, 1 chamber 17ml volume		X		X	
5316-1085	Uno sampler, V4A/PTFE, 85cm, 1 chamber 17ml volume		X		X	
5316-1150	Uno sampler, V4A/PTFE, 150cm, 1 chamber 17 ml vol.		X		X	
5316-2055	All-layer-sampler, V4A, 55 cm, 3 openings, 70 ml		X			
5316-2085	All-layer-sampler, V4A, 85 cm, 3 openings, 130 ml		X			
5316-2150	All-layer-sampler, V4A, 150 cm, 5 openings, 260ml		X			
5316-3055	Multi sampler, V4A, 55cm, 3 chambers per 14ml vol.		X			
5316-3085	Multi sampler, V4A, 85cm, 3 chambers per 17ml vol.		X			
5316-3150	Multi sampler, V4A, 150cm, 5 chambers per 17 ml		X			
5317-0030	Emptying device, V2A, length 100 cm, complete	X				
5317-0033	Emptying funnel separate, V2A	X				
5317-1058	Multi sampler, Alu, 55 cm, 3 chambers per 14 ml			X		
5317-1085	Multi sampler, Alu, 85 cm, 3 chambers per 17 ml			X		
5317-1150	Multi sampler, Alu, 150cm, 5 chambers per 17 ml			X		
5317-2055	All-layer-sampler, Alu, 55 cm, 3 openings, 70 ml			X		
5317-2085	All-layer-sampler, Alu, 85 cm, 3 openings, 130 ml			X		
5317-2150	All-layer-sampler, Alu, 150 cm, 5 openings, 260ml			X		
5318-0145	Zone sampler Spiralus, Alu/PA, 145cm, 770ml volume			X		X
5318-0200	Zone sampler Spiralus, Alu/PA, 200cm, 1100 ml vol.			X		X
5319-1085	Jumbo sampler, all layer, Alu, 85cm, 880 ml volume			X		
5319-1150	Jumbo sampler, all layer, Alu, 150 cm, 1700ml vol.			X		
5319-1250	Jumbo sampler, all layer, Alu, 250 cm, 2900ml vol.			X		
5319-2085	Jumbo sampler, closed chambers, Alu/PTFE, 85 cm			X	X	

Article	Description	V2A	V4A	Alu	PTFE	PA
5319-2150	Jumbo sampler, closed chambers, Alu/PTFE, 150 cm			X	X	
5319-2250	Jumbo sampler, closed chambers, Alu/PTFE, 250 cm			X	X	
5320-0150	Sample-spoon, V2A, 150 mm, 1,5 ml, autoclavable	X				
5320-0180	Sample-spoon, V2A, 180 mm, 2 ml, autoclavable	X				
5320-0210	Sample-spoon, V2A, 210 mm, 6 ml, autoclavable	X				
5320-0300	Sample-spoon, V2A, 300 mm, 6 ml, autoclavable	X				
5320-0400	Sample-spoon, V2A, 400 mm, 6 ml, autoclavable	X				
5320-0500	Sample-spoon, stainless steel, V2A, 500 mm, 9 ml	X				
5321-0300	Sampling lance Milky, V4A, 480 mm, 150 ml, Ø 28 mm		X			
5322-1300	QuickPicker V4A, 300 mm, 75 ml volume, ISTA-Norm		X			
5323-0200	Ice borer, V4A, inner-xouter-Ø 16x21mm, depth 20cm		X			
5323-2010	Ice borer, V4A, inner-xouter-Ø 50x54mm, depth 20cm		X			
5324-0001	PharmaSpoon, V4A, length 650 mm, hand polished		X			
5324-0010	PharmaScoop, V4A, 10 ml vol., Ø 25 mm, length 45mm		X			
5324-0011	PharmaScoop, V4A, 50 ml vol., Ø 38 mm, length 70mm		X			
5324-0012	PharmaScoop, V4A, 100ml vol., Ø 50mm, length 100mm		X			
5324-0013	PharmaScoop, V4A, 750ml vol., Ø 95mm, length 210mm		X			
5324-0014	PharmaScoop, V4A, 1000 ml vol., Ø110 mm, 235 mm		X			
5324-0015	PharmaScoop, V4A, 1250 ml vol., Ø125 mm, 195 mm		X			
5324-0016	PharmaScoop, V4A, 1700 ml vol., Ø135 mm, 230 mm		X			
5324-0017	PharmaScoop, V4A, 2500 ml vol., Ø185 mm, 215 mm		X			
5324-0018	PharmaScoop, V4A, 500ml vol., Ø 95mm, length 140mm		X			
5324-0019	PharmaScoop with long handle, V4A, 1000 ml		X			
5326-0001	Liquid-Sampler, open with thumb, V4A, 50 ml		X			
5326-0002	Liquid-Sampler, open with thumb, V4A, 100 ml		X			
5326-0003	Liquid-Sampler, place onto bottle, V4A, 50 ml		X			
5326-0004	Liquid-Sampler, place onto bottle, V4A, 100 ml		X			
5334-8000	Stainless steel measuring spoons set (1,25-15 ml)	X				
5338-0410	StickProof, V4A, length 410 mm, Ø 25mm, ISTA-Norm		X			
5338-1100	StickProof, V4A, length 1100mm, Ø 25mm, ISTA-Norm		X			
5339-0300	IceSampler, V4A, length 300mm, vol. 50ml, Ø 38mm		X			
5339-0550	IceSampler, V4A, length 550mm, vol. 50ml, Ø 38mm		X			
5339-1050	IceSampler, V4A, length 1050mm, vol. 50ml, Ø 38mm		X			
5340-1100	SiloPicker, V4A, vol. 100 ml, 130 cm, outer-Ø 50mm		X			
5341-1002	Suction head, suction drill 2mm, Mini ViscoSampler		X			
5341-1004	Suction head, suction drill 4mm, Mini ViscoSampler		X			
5341-1006	Suction head, suction drill 6mm, Mini ViscoSampler		X			
5341-1060	Mini ViscoSampler, V4A, 60 cm, vol. 30 ml, Ø 15 mm		X			
5341-1100	Mini ViscoSampler, V4A, 100 cm, vol. 50 ml, Ø 15mm		X			
5341-1150	Mini ViscoSampler, V4A, 150 cm, vol. 75 ml, Ø 15mm		X			
5346-1000	PharmaPicker, V4A, w/collection cylinder 0,1-1,2ml		X			
5346-1001	Collection cylinder for 0,1-1,2ml, w/o head		X			
5346-1010	Volume head, capacity 0,1 ml, PharmaPicker		X			
5346-1020	Volume head, capacity 0,2 ml, PharmaPicker		X			
5346-1060	Volume head, capacity 0,6 ml, PharmaPicker		X			
5346-1100	Volume head, capacity 1 ml, PharmaPicker		X			
5346-2000	PharmaPicker, V4A, w/coll. cylinder for 1,25-2,5ml		X			
5346-2001	Collection cylinder for 1,25-2,5ml, w/o head		X			
5346-2200	Volume head, capacity 2 ml, PharmaPicker		X			

Article	Description	V2A	V4A	Alu	PTFE	PA
5346-2250	Volume head, capacity 2,5 ml, PharmaPicker		X			
5347-0005	Scratcher, V4A, 5 ml, flute length 22 mm, Ø 25 mm		X			
5347-0010	Scratcher, V4A, 10 ml, flute length 45 mm, Ø 25 mm		X			
5347-0025	Scratcher, V4A, 25 ml, flute length 110 mm, Ø 25mm		X			
5351-1000	All-layer-scratcher, V4A, Ø 25mm, 100cm, vol.180ml		X		X	
5354-1000	Stainless steel scoop, V2A, 115cm, capacity 1000ml	X				
5354-2000	Stainless steel beaker, V2A, 1000 ml, TeleScoop	X				
5354-3000	Pendulum beaker V2A, capacity 1000 ml, TeleScoop	X				
5356-1500	EasyScoop, aluminium, 1500 ml, inner grip			X		
5356-2000	EasyScoop, aluminium, 2000 ml, inner grip			X		
5357-0001	Sampler Liquid CupSampler, V4A, 100 ml, 100 cm		X			
5358-0000	Sample insert, 0,0 ml, Novartis Multi		X			
5358-0005	Sample insert, 0,5 ml, Novartis Multi		X			
5358-0010	Sample insert, 1 ml, Novartis Multi		X			
5358-0015	Sample insert, 1,5 ml, Novartis Multi		X			
5358-0020	Sample insert, 2 ml, Novartis Multi		X			
5358-0100	Sample insert, 0,0 ml, Novartis Uno		X			
5358-0105	Sample insert, 0,5 ml, Novartis Uno		X			
5358-0110	Sample insert, 1 ml, Novartis Uno		X			
5358-0120	Sample insert, 2 ml, Novartis Uno		X			
5358-0130	Sample insert, 3 ml, Novartis Uno		X			
5358-0140	Sample insert, 4 ml, Novartis Uno		X			
5358-0600	Novartis Uno, V4A, 600 mm, 2 chambers, w/o inserts		X			
5358-1000	Novartis Multi, V4A, 100cm, 14 chambers, w/o ins.		X			
5358-1001	Novartis Uno, V4A, 1000mm, 2 chambers, w/o inserts		X			
5358-1400	Novartis Uno, V4A, 1400mm, 2 chambers, w/o inserts		X			
5359-1600	Valve lance, V4A, length 550 mm, Ø 32 mm		X			
5359-1800	Valve lance, V4A, length 800 mm, Ø 32 mm		X			
5360-0600	EasySampler V4A 60ml		X			
5360-1000	EasySampler, V4A, length 104,5 cm, volume 200 ml		X			
5369-0300	Pigment lance, flute open.25%, length 300mm, 145ml		X			
5369-0310	Pigment lance, flute open. 50%, flute 300mm, 85 ml		X			
5369-0600	Pigment lance, flute open. 25%, flute 600mm, 312ml		X			
5369-0610	Pigment lance, flute open. 50%, flute 600mm, 179ml		X			
5369-0900	Pigment lance, flute open. 25%, flute 900mm, 480ml		X			
5369-0910	Pigment lance, flute open. 50%, flute 900mm, 273ml		X			
5370-0120	Scoop AISI 304, V2A, capacity 120 ml, Ø 50 mm	X				
5370-0450	Scoop AISI 304, V2A, capacity 450 ml, Ø 75 mm	X				
5370-1000	Scoop AISI 304, V2A, capacity 1000 ml, Ø 100 mm	X				
5370-1500	Scoop AISI 304, V2A, capacity 1500 ml, Ø 120 mm	X				
5370-2000	Scoop AISI 304, V2A, capacity 2500 ml, Ø 150 mm	X				
5371-0100	FoodScoop, V2A, 100ml, total l:215mm, LxW 120x60mm	X				
5371-0200	FoodScoop, V2A, 200ml, total l:230mm, LxW 140x75mm	X				
5371-0250	FoodScoop, V2A, 250ml, total l:245mm, LxW 155x85mm	X				
5371-0350	FoodScoop, V2A, 350ml, total: 290mm, LxW 180x100mm	X				
5371-0500	FoodScoop, V2A, 500ml, total: 350mm, LxW 220x130mm	X				
5380-0508	SiloPicker extension 50 cm		X			
5380-1008	SiloPicker extension 100 cm		X			
5382-0040	Sampler Tubus, V4A, 40 cm, Ø 40 mm, ISTA standard		X			

Article	Description	V2A	V4A	Alu	PTFE	PA
5382-0080	Sampler Tubus, V4A, 80 cm, Ø 40 mm, ISTA standard		X			
5383-1000	Spatula stainless steel, V4A, 300 mm		X			
5386-0001	Micro spatula stainless steel, V2A, length 170 mm	X				
5386-0010	Palette knife spatula stainl. steel, LxW 200x21 mm	X				
5386-0012	Palette knife spatula stainl. steel, LxW 250x25 mm	X				
5386-0022	Spoon spatula stainless steel, length 180 mm, 4 ml	X				
5386-0025	Spoon spatula stainless steel, length 300 mm, 11ml	X				
5386-0100	Forceps, stainless steel, sharp, bent form, 105 mm	X				
5386-0102	Forceps, stainless steel, sharp, bent form, 130 mm	X				
5386-0104	Forceps, stainless steel, sharp, bent form, 160 mm	X				
5386-0200	Forceps, stainless steel, sharp, straight, 105 mm	X				
5386-0202	Forceps, stainless steel, sharp, straight, 130 mm	X				
5386-0204	Forceps, stainless steel, sharp, straight, 160 mm	X				
5386-0300	Forceps, stainless steel, blunt, straight, 105 mm	X				
5386-0302	Forceps, stainless steel, blunt, straight, 130 mm	X				
5386-0304	Forceps, stainless steel, blunt, straight, 160 mm	X				
5388-1000	BeefSteaker, V2A, Ø 25mm	X				
9604-0001	Funnel, stainless steel, V2A, outer-Ø 60 mm	X				
9604-0003	Funnel, stainless steel, V2A, outer-Ø 100 mm	X				
9604-0005	Funnel, stainless steel, V2A, outer-Ø 150 mm	X				
9604-0008	Trichter Edelstahl, V2A, Außen-Ø 300 mm	X				