ZERTIFIKAT
CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft
Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science

6202-xxxx

All-purpose tube, w/ screw cap, PP, transp.,various size

Declaration of Conformity

Food Contact Materials

General Compliance

We hereby confirm that the above-mentioned product complies with the German Commodity Regulation, with the Regulation (EU) No. 10/2011 and subsequent amendments as well as the EU-Framework Regulation (EC) No. 1935/2004.

Conditions of use / Compliance with threshold values

a) Specifications of the intended use or limitations:
- Type(s) of food intended to come into contact with the product/material: aqueous, acidic, fatty, alcohol (up to 10% alcohol) and solid food
- Duration and temperature of treatment and / or storage while in contact with the food: Long term contact / - storage 10d / 40°C
- Ratio of the area of the food contact material to the volume used to determine the compliance of the plastic food contact material or article: 6 dm2/kg
- Test requirements: Food stimulants A, B, D2 (10% Ethanol; 3% Acetic Acid; Oil)

b) The global migration values are lower than the admitted threshold values. The tests are performed according to Article 17/18 of Regulation (EU) No. 10/2011 in conjunction with Annex V.
Information about restricted substances and or specifications / compliance of the limits is confirmed.

Stoffbezeichnung / Name	CAS-Nr.	Beschränkung / Limit [mg/kg]
1,2,3-Trideoxy-4,6:5,7-bis-O- ((4-propylphenyl)methylene)-nonitol	882073-43-0	5
Acetonitril	7128-64-5	0,6
N,N-Bis(2-hydroxyethyl)alkyl(C8-C18)amin	70955-14-5	1,2
Aluminium		1,0

Information about dual use additives according to Regulation (EU) No. 10/2011

Stoffbezeichnung / Name	CAS-Nr.
Glycerinmonostearat (Glycerol Monostearate)	31566-31-1 / E 471
Benzoesäure (Benzoic Acid)	CAS No.65-85-0 / Ref.-Nr. 37600 / E 210

No usage of a plastic functional barrier.

ZERTIFIKAT
CERTIFICATE

Pumpen | Probenehmer | Laborbedarf
für Labor, Industrie und Wissenschaft

Pumps | Sampling | Plastic Labware
for Laboratory, Industry, Science

6202-xxxx

All-purpose tube, w/ screw cap, PP, transp., various size

In terms of Regulations (EU) No. 10/2011, (EC) No. 1935/2004 and (EC) No. 2023/2006 there are no objections against the use of the product for the manufacture of articles intended for food contact. The raw materials used comply with the demands of the Regulation (EU) No. 10/2011. This Regulation provides guidelines for the selection of test conditions to be used for various food products. This declaration is valid for the product delivered by us as specified above. According to that and under consideration of the food contact conditions stated, the product complies with the stipulations of these Regulations regarding the packaging of food products to be packed. The user shall verify himself that the product is suitable for the intended food to be packed beyond the stipulations of the Directives. In particular it is emphasized that in case of printed material, no contact is allowed between the printing ink and the food product.

Further tests are not carried out by us.

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of certificate is 002284-QM15.

September 21, 2026

Bürkle GmbH, Bad Bellingen
Martin Saint-Denis, Managing Director